

## Dry & Wet Sieve Shaker Machine

A small and light-weight laboratory separator capable of separating up to 8 sieves. Although small in size, this shaker still utilizes three dimensional vibration for maximum performance and complete separation



### Specification:

|                  |                  |              |
|------------------|------------------|--------------|
| Model            | EI-S200          | EI-S300      |
| Number of Sieves | 8                |              |
| Sieve Diameter   | 200mm            | 300mm        |
| Feed Capacity    | 200g             |              |
| Timer            | 0 ~ 99min59s     |              |
| Maximum Speed    | 1420rpm/min      |              |
| Noise            | 60db             |              |
| Power            | 220V, 50Hz       |              |
| Dimension        | 360 x 300 x736mm |              |
| Weight           | Approx. 25kg     | Approx. 31kg |