

Metallographic Specimen grinding & Polishing Machine



Application

Metallographic Specimen grinding Machine is suitable for the rough grinding, and fine & half fine grinding of the metallographic specimen

Equipped with water inlet which can realize the cooling onto the specimen to prevent damaged caused by the heat.

This machine features convenient operation, stable and reliable performance, it's ideal sample preparing equipment to be used in factories, R&D and university laboratories.

Specification

Brand	ECO
Model	EI-G1
Disc Diameter	200 mm
Rotation Speed	900 rpm
Abrasive paper Diameter	200 mm
Rate Voltage	220-240V, 50Hz
Dimension	400 x 470 x 550 mm
Weight	Approx. 15kg