

Manual Sample Press Machine



This manual sample press machine is to punch a test specimen from a formed sheet, such as plastic (Elastomer), film, vulcanized rubber, paper, fabrics, etc. suitable for various type of dumbbell, circle, rectangle shape of die blocks, user friendly, easily operate this manual press by one hand.

Specification:

Model	EI-P20M
Cutting Area	100cm ²
Maximum Stroke	20mm
Vertical Space	45mm
Platen	Ø120 x 10mm
Working Table	Ø135mm
Dimension	260 x 170 x 420mm
Weight	Approx. 8kg
Accessories	Tools x 1 Cutting Board x 1 Compression Platen Die Block - 100cm ²