

# Notch Impact Sample Cutter



## Descriptions:

A semi-auto preparation notched specimens of the non-metallic material impact toughness test.

## Features:

- Fast and accurate with multi-sample loading
- Motorized operation
- User Friendly & safety guard

## Specifications:

Model	EI-NM
Notched	A type: 45°, r = 0.25mm B type: 45°, r = 1.0mm C type: 45°, r = 0.1mm
Travel Length	160mm
Cutting Speed	150rpm/min
Standard Accordance	ISO179, ISO180, ASTM D256
Power	240V, 50Hz
Dimension	600 x 270 x 300mm
Weight	Approx. 30kg